

- In addition to routine checks for each use, PPE should regularly undergo a detailed inspection by a competent person. Petzl recommends an inspection every 12 months and after any exceptional event in the life of the product.
 - PPE inspection should be conducted with the manufacturer's Instructions for Use.
- Download the instructions at [PETZL.COM](https://www.petzl.com).

NEST - STEF

User information

Surname

Address

Identifier

PPE information

Model

Serial number

Year of manufacture

Date of purchase

Date of first use

Manufacturer: Petzl, ZI Cidex 105A - 38920 Crolles - France

 Good condition (G)	 To monitor (TM)	 To repair (TR)	 Do not use, retire (R)	N/A Not applicable
---	--	---	---	-----------------------

1. Known product history

Usage conditions or exceptional event during use (examples: fall or fall arrest, use or storage at extreme temperatures, modification outside manufacturer's facilities...):

				N/A
---	---	---	---	-----

2. Preliminary observations

- Verify the presence and legibility of the serial number and the CE mark.
- Verify that the product lifetime has not been exceeded.
- Compare with a new device to verify there are no modifications or missing elements.

3. Checking the condition of the STEF

- Check the condition of the connectors: see the inspection form for your connector model at [PETZL.COM](https://www.petzl.com)
- Check the condition of the swivel: see the inspection form for Rigging equipment at [PETZL.COM](https://www.petzl.com)
- Check the condition of the straps (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals).
- Check the condition of the safety stitching on both sides. Look for any threads that are loose, worn, or cut. Safety stitching is identified by thread of a different color than that of the webbing.
- Check the condition of the metal ring (deformation, cracks, marks, wear, corrosion).

NEST

4. Checking the external elements

- Check the condition of the carry handles (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals).
- Check the condition of the flap closure straps (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals) and the function of the fastening buckles.
- Check the condition of the flaps, the head and foot covers and their adjustment straps (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals) and the function of the fastening buckles.
- Check the condition of the tarp's exterior (deformation, cracks, marks, wear).
- Check the condition of the rivets (loose or missing rivet, deformation, cracks, marks, wear, corrosion).
- Remove and check the condition of the longitudinal stays (deformation, cracks, marks, wear, corrosion).
- Check the condition of the accessory attachment points (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals).
- Remove and check the condition of the transverse stays (deformation, cracks, marks, wear, corrosion).

5. Checking the attachment points

- Check the condition of the metal attachment points (deformation, cracks, marks, wear, corrosion).
- Check the condition of the straps and the textile attachment points (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals).
- Check the condition of the safety stitching on both sides. Look for any threads that are loose, worn, or cut. Safety stitching is identified by thread of a different color than that of the webbing.

6. Checking the condition of the patient positioning straps.

- Check the condition of the straps (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals).
- If necessary, lift the foam pad.
- Verify that items are present on the strap ends.
- Check the condition of the safety stitching on both sides. Look for any threads that are loose, worn, or cut. Safety stitching is identified by thread of a different color than that of the webbing.
- Check the condition of the FAST adjustment buckles (deformation, cracks, marks, wear, corrosion) and that they are working properly.
- Check the condition of the foot supports (cuts, swelling, damage and wear due to use, to heat, to contact with chemicals) and that the adjustment buckles are working properly.
- Check the condition of the foam pad.

Comments (detail here any defects found on the product and accomplished tasks)

VERDICT

Product **fit** to remain in service

Product **unfit** to remain in service

Inspected by

Company

Date

Next inspection date